

Removing and Replacing the Transducer Battery



Correct use




Incorrect use



Recycle according to local laws

CAUTION: NEVER immerse a transducer in liquid if the battery drawer is open or has been removed. Check the condition of the drawer before refitting, and verify that it is fully closed before use.



8




9



10



Wipe off excess silicone grease.



Informaton for Silicone Grease

Manufacturer:

ELBE

L. Böwing GmbH

Delkenmheimer Str. 18

65719 Hofheim a. Ts.

Germany

Product: Elbesil BM Silicone Grease

Lot number: Impressed on the fold at the bottom of the tube

Minimum shelf life: See sticker on battery exchange kit packaging

